


По вопросам продаж и поддержки обращайтесь:

Архангельск (8182)63-90-72	Краснодар (861)203-40-90	Рязань (4912)46-61-64
Астана (7172)727-132	Красноярск (391)204-63-61	Самара (846)206-03-16
Белгород (4722)40-23-64	Курск (4712)77-13-04	Санкт-Петербург (812)309-46-40
Брянск (4832)59-03-52	Липецк (4742)52-20-81	Саратов (845)249-38-78
Владивосток (423)249-28-31	Магнитогорск (3519)55-03-13	Смоленск (4812)29-41-54
Волгоград (844)278-03-48	Москва (495)268-04-70	Сочи (862)225-72-31
Вологда (8172)26-41-59	Мурманск (8152)59-64-93	Ставрополь (8652)20-65-13
Воронеж (473)204-51-73	Набережные Челны (8552)20-53-41	Тверь (4822)63-31-35
Екатеринбург (343)384-55-89	Нижний Новгород (831)429-08-12	Томск (3822)98-41-53
Иваново (4932)77-34-06	Новокузнецк (3843)20-46-81	Тула (4872)74-02-29
Ижевск (3412)26-03-58	Новосибирск (383)227-86-73	Тюмень (3452)66-21-18
Казань (843)206-01-48	Орел (4862)44-53-42	Ульяновск (8422)24-23-59
Калининград (4012)72-03-81	Оренбург (3532)37-68-04	Уфа (347)229-48-12
Калуга (4842)92-23-67	Пенза (8412)22-31-16	Челябинск (351)202-03-61
Кемерово (3842)65-04-62	Пермь (342)205-81-47	Череповец (8202)49-02-64
Киров (8332)68-02-04	Ростов-на-Дону (863)308-18-15	Ярославль (4852)69-52-93

Единый адрес: pca@nt-rt.ru **Веб-сайт:** www.precisa.nt-rt.ru

Весы 320 XB Gold/Carat PRECISA



SWISS MADE 

Gold/Carat Balances Series 320 XB

Precisa gold and carat balances are precision objects as valuable as gold

Precisa gold and carat balances owe its extremely high reliability and durability to the well-known competence of Precisa's engineering and, as a result, exceed all established standards. These balances were developed for professional use in a wide range of jewelry and precision weighing applications.



- Top Swiss quality
- Elegant compact design
- Practical and easy to use
- Quick response and high stability
- Wide range of models
- High-contrast vacuum fluorescent display
- High resolution balance for an attractive price
- 3 doors draft shield for the 0.001 ct balance
- SCS fully automatic self-calibration system (optional)
(time and temperature controlled internal calibration)
- RS232 interface for PC and Printer
- Two-years manufacturer's warranty




The value of a gem is not only in its look!

Elegance and beauty in the design are the unmistakable external features. Precision, solid quality and problem-free operation is the essential internal merit. That is why in every respect Precisa gold and carat balances are precision objects as valuable as gold. Except for the price!

Precisa

■ The Balance of Quality ■



5 keys for uncomplicated operation

The requirement for perfect functioning is easy operation. Clear user guidance guarantees error-free weighing.

Exclusive: Precisa anti-theft device

A gem like this balance must be protected.

Precisa balance is fitted with the option of a mechanical and, exclusively, an electronic anti-theft device which automatically deactivates the balance. It can be switched on again using the secret personal code known only to you.

Precisa Gold/Carat Balances Series 320: Technical Specification

Model	Capacity	Readability	Repeatability	Linearity	Pan Size (mm/inch)
XB 600M-C *	600 ct / 120 g	0.001 ct / 0.001 g	0.001 ct / 0.001 g	0.001 ct / 0.001 g	Ø 80 / 3.1
XB 1100M-C *	1 100 ct / 220 g	0.001 ct / 0.001 g	0.001 ct / 0.001 g	0.001 ct / 0.001 g	Ø 80 / 3.1
XB 3100C-C *	3 100 ct / 620 g	0.01 ct / 0.01 g	0.01 ct / 0.01 g	0.01 ct / 0.01 g	135 x 135
XB 620C-G *	620 g	0.01 g	0.01 g	0.01 g	170 x 170
XB 1200C-G *	1 200 g	0.01 g	0.01 g	0.01 g	170 x 170
XB 2200C-G *	2 200 g	0.01 g	0.01 g	0.01 g	170 x 170
XB 4200C-G *	4 200 g	0.01 g	0.01 g	0.01 g	170 x 170

* models available with EC/OIML type approval

Additional Application & Features	Accessories	Series 320
Further special accessories and options on demand	Diamond weighing pan	350-8322

Features and Applications



Swiss made



Support IQ/OQ/PQ



Density determination



Animal Weighing



Vacuum Fluorescent Display



Anti theft protection password and mechanical



Counting program



Automatic calibration warning



% weighing



Units g, mg, oz, lb, tael, etc.



Interface RS232 for PC printer



Below balance weighing



Net total application
Add-up summation
Dosage program

По вопросам продаж и поддержки обращайтесь:

Архангельск (8182)63-90-72	Краснодар (861)203-40-90	Рязань (4912)46-61-64
Астана (7172)727-132	Красноярск (391)204-63-61	Самара (846)206-03-16
Белгород (4722)40-23-64	Курск (4712)77-13-04	Санкт-Петербург (812)309-46-40
Брянск (4832)59-03-52	Липецк (4742)52-20-81	Саратов (845)249-38-78
Владивосток (423)249-28-31	Магнитогорск (3519)55-03-13	Смоленск (4812)29-41-54
Волгоград (844)278-03-48	Москва (495)268-04-70	Сочи (862)225-72-31
Вологда (8172)26-41-59	Мурманск (8152)59-64-93	Ставрополь (8652)20-65-13
Воронеж (473)204-51-73	Набережные Челны (8552)20-53-41	Тверь (4822)63-31-35
Екатеринбург (343)384-55-89	Нижний Новгород (831)429-08-12	Томск (3822)98-41-53
Иваново (4932)77-34-06	Новокузнецк (3843)20-46-81	Тула (4872)74-02-29
Ижевск (3412)26-03-58	Новосибирск (383)227-86-73	Тюмень (3452)66-21-18
Казань (843)206-01-48	Орел (4862)44-53-42	Ульяновск (8422)24-23-59
Калининград (4012)72-03-81	Оренбург (3532)37-68-04	Уфа (347)229-48-12
Калуга (4842)92-23-67	Пенза (8412)22-31-16	Челябинск (351)202-03-61
Кемерово (3842)65-04-62	Пермь (342)205-81-47	Череповец (8202)49-02-64
Киров (8332)68-02-04	Ростов-на-Дону (863)308-18-15	Ярославль (4852)69-52-93

Единый адрес: pca@nt-rt.ru **Веб-сайт:** www.precisa.nt-rt.ru